


Down-Conversion and Photon Counting

Tristan Perry & Jaden Miller


Potential Application

- Optimizing Photon counts
 - Speed of light through some medium
 - Laser testing
 - Detection of downconversion
- Testing downconversion
 - Entanglement experiments
 - Probability testing of photon location



Electromagnetism and the Photon

- Light as a wave and a particle
- Classical Period
 - Newton's "corpuscles"
 - Thomas Young's Double Slit
 - Malus's Law
 - Maxwell
- Modern Period
 - Planck's Quantum Energy & Frequency, $E = hf$
 - Einstein's Photoelectric theory
 - 2022 Nobel Prize






EM wave propagation:

- $\mathbf{E} = E_0 e^{i(\omega t - \mathbf{k} \cdot \mathbf{r})}$
- $\mathbf{H} = H_0 e^{i(\omega t - \mathbf{k} \cdot \mathbf{r})}$

https://em.geosci.xyz/content/maxwell1_fundamentals/harmonic_planewaves_homogeneous/impedancephase.html

Photons as a Quantum Particle

- *A hands-on introduction to single photons and quantum mechanics for undergraduates*, Dickinson College PA
 - Brett J. Pearson and David P. Jackson




- *Hologram of a single photon*, Warsaw University
 - Radosław Chrapkiewicz, Michał Jachura, Konrad Banaszek & Wojciech Wasilewski

[\[1509.02890\] Hologram of a Single Photon \(arxiv.org\)](#)

Counting “Single” Down-converted Photons

Objective: Optimizing the count of single photons



Theory: Detecting Down-converted Photons

$$E = hf = hc/\lambda$$

$$E_1 = E/2 = hf/2 = hc/2\lambda$$

$$E_2 = E/2 = hf/2 = hc/2\lambda$$


Desired λ to be detected is twice that of the StellarPro:

$$\lambda_{SP} = (514.5 \pm 1\%) \text{nm}$$


[Stellar-Pro ML-150.pdf \(modu-laser.com\)](#)

$$\lambda_{DC} = 2(514.5 \pm 1\%) \text{nm} = (1029 \pm 2\%) \text{nm}$$


Conservation of Energy



Conservation of Momentum




Theory: BBO Downconversion



Setup & Alignment

- ModuLaser StellarPro ML/150 Ar laser ($\lambda=514.5\text{nm}$)
- UniPhase Helium Neon Gas Laser 1508P-0
- EKSPA Optics BBO 6x6x1 th29.2 ph0 (Beta Barium Borate)
- Photon Detectors (IDQ-id100)



Setup & Procedure


- Oscilloscope (Tektronix TDS 2024C)
- TeachSpin Pulse Counters/Interval Timers
- ThorLabs waveplate and iris
- IDQ-id100/150 UN single Photon detectors

$$V_{\text{out}} = (3.3 \pm 1\%)V \text{ per pulse}$$



Procedure


1. Filters to prevent overload
2. Count data taken with and without BBO crystals
3. 30 trials to give avg



Data & Results


- For $N = 30$

Detector	W_{AVG} Photon count with BBO	δW_{AVG}
Left	2318.7	8.7
Right	2059.3	8.2
Both	2138.3	6.0
	W_{AVG} Photon count No BBO	
Left	48.7	1.3
Right	37.7	1.1




Errors, Improvements

1. Laser Variance
2. False Counts from Signal Noise - 20 photons background S:N
 - a. Equipment & Environment
3. Visual Setup Alignment / Detector Alignment / Realignment
4. Detector limitations




Probability versus λ



Improvements

Future experiments include:

1. IDQ IQ230 IR detector (900-1700nm)
2. Potential entanglement pair detection



References

Hossieni, H., Fatah, J. M. A., Mohammad, S., & Naby, M. (2018). From Hyperion to photon, a brief survey in the timeline of photon. *Physics Essays*, 31(3), 246–253. <https://doi.org/10.4006/0836-1398-31.3.246>

Geometry Matters. (2022, April 12). *Hologram of a single photon*. Geometry Matters. Retrieved November 30, 2022, from <https://geometrymatters.com/hologram-of-a-single-photon/>

Pearson, B. J., & Jackson, D. P. (2010, February 10). *A hands-on introduction to single photons and quantum mechanics for undergraduates*. <https://www.aapt.org/>. Retrieved November 20, 2022.

Dehlinger, D., & Mitchell, M. W. (1970, January 1). *Entangled photon apparatus for the Undergraduate Laboratory*. American Association of Physics Teachers. Retrieved December 2, 2022, from <https://aapt.scitation.org/doi/10.1119/1.1498859>

Galvez, K. (2022). *Alpha immersion – hands-on Activities*. Colgate University. Retrieved December 2, 2022, from <http://egalvez.colgate.domains/pql/wp-content/uploads/2019/08/Activities19.pdf>


Questions

Special Thanks to:
Dr. Kiko Galvez
Dr. Brooke Hester
Dr. Patricia Allen
App State IT department
TeachSpin help desk
Jaden Miller

Additional info: Photodiodes

- Photons incident on reverse biased PN junction
- Photons break bonds to create e^- - hole pairs
- $hf > E_{bg}$
- Voltage produced relative to intensity

$$(h = 6.62607015 \times 10^{-34} \text{ m}^2 \text{ kg} / \text{s})$$



Additional Info: Video of Photon from MIT

https://www.youtube.com/watch?v=_QAPQO6EL8o